

Калиев йодид

extra pure, SLR

Промо цена!
119 лв/кг

Разфасовка
5 кг



Specification Sheet

Product: Potassium iodide

Code: P/5840

Grade: SLR

Application: Specified reagent for general lab work

Expiry Phrase: Use within 5 years of opening

Test	Guaranteed	Units
Assay	≥ 99.5	%
Calcium (Ca)	≤ 20	ppm
Copper (Cu)	≤ 10	ppm
Iodate (IO ₃)	≤ 10	ppm
Iron (Fe)	≤ 10	ppm
Lead (Pb)	≤ 10	ppm
Loss on drying	≤ 0.5	%
Magnesium (Mg)	≤ 20	ppm
Sodium (Na)	≤ 2000	ppm
Total phosphorus (P)	≤ 50	ppm
Total silicon (Si)	≤ 50	ppm
Total sulfur (S)	≤ 100	ppm
Zinc (Zn)	≤ 50	ppm

Additional Information:

Mrs K S Cluskey BSc
Quality & Compliance Manager.

Thursday 03rd of January 2013 11:44:03 AM
Any queries please phone: +44(0)1509 555158.

ThermoFisher
SCIENTIFIC

The world leader in serving science

 **Fisher Chemical**

Лабимекс ООД

Адрес: София 1404, ж.к. „Гоце Делчев“ бл. 261, офис "Ч"

Телефони: 02 422 41 60

Тел./ факс: 02 962 15 01

E-mail: office@labimex.bg

www.labimex.bg